YNK MEDIKAL - HEALING ABUTMENT INSTRUCTIONS FOR USE (IFU)

ISO 13485:2016 & CE Certified Dental Implant Prosthetic Component Guide

1. INTENDED USE / KULLANIM AMACI:
   Designed for peri-implant soft tissue contouring during the transmucosal healing phase.
2. MATERIAL / HAMMADDE:
   Medical Grade 5 Titanium Alloy (Ti-6Al-4V ELI - ASTM F136 / ISO 5832-3).
3. RECOMMENDED TORQUE / ONERILEN SIKMA TORKU:
   10 - 15 Ncm (Hand-tighten or calibrated prosthetic torque wrench with 1.20/1.25 mm hex driver).
4. STERILIZATION PROTOCOL / STERILIZASYON PROTOKOLU:
   Components are supplied NON-STERILE. Prior to intraoral placement, clean ultrasonically
   and steam-autoclave in a medical pouch at 134 deg C for a minimum of 5 minutes.
5. SINGLE PATIENT USE / TEK HASTA KULLANIMI:
   Do not reuse across multiple patients due to biological cross-contamination and thread wear.

Manufacturer: YnK Medikal Mekanik Elektronik San. ve Tic. Ltd. Sti.
Ikitelli OSB Mah. Marmara Sanayi Sitesi R-Blok 10/B Istanbul / Turkiye | info@ynkmedikal.com