YNK MEDIKAL - MULTI-UNIT, BALL ATTACHMENT & LOCATOR GUIDE

Full-Arch Screw-Retained & Implant-Supported Overdenture Prosthetic Protocol

1. MULTI-UNIT ABUTMENT SYSTEM:
   Engineered for passive-fit screw-retained bridges and hybrid full-arch restorations.
   - Fixture-level Multi-Unit Seating Torque: 30 Ncm (Mini platforms: 25 Ncm).
   - Prosthetic Coping Screw (M1.4 / M1.6 / M1.8): 15 Ncm.
2. BALL ATTACHMENT (D: 2.50 mm) & LOCATOR OVERDENTURE ABUTMENTS:
   Available in gingival heights GH 1.0 mm, 2.0 mm, 3.0 mm, 4.0 mm, and 5.0 mm.
   Select gingival cuff height matching mucosal thickness + 0.5 mm supragingival emergence.
   Recommended seating torque: 25 - 30 Ncm.
3. STERILIZATION: Steam autoclave at 134 deg C for 5 minutes before clinical delivery.

Manufacturer: YnK Medikal Mekanik Elektronik San. ve Tic. Ltd. Sti. | +90 212 494 37 36