YNK MEDIKAL - CEMENTABLE & ANGLED ABUTMENT TORQUE PROTO

Straight, 15 deg & 25 deg Titanium Abutment Clinical & Laboratory Guide

1. COMPATIBLE PLATFORMS / UYUMLU BAGLANTI SISTEMLERI:
   Cross Standard & Mini, Octagon, Internal Hex 2.42mm, Conical Hex Standard & Mini.
2. RECOMMENDED CLINICAL TORQUE / KLINIK VIDA TORK DEGERLERI:
   - Cross Standard / Internal Hex / Conical Hex Standard Abutment Screw: 30 Ncm
   - Cross Mini / Conical Hex Mini Narrow Platform Screw: 20 - 25 Ncm
   - Prosthetic Retaining Screw (Multi-Unit Cap Screw): 15 Ncm
3. PREPARATION & STERILIZATION:
   Perform any extraoral supragingival milling with copious water irrigation and carbide/titanium burs. Steam sterilize at 134 deg C for 5 minutes prior to final intraoral seating.

Manufacturer: YnK Medikal Mekanik Elektronik San. ve Tic. Ltd. Sti. | ynkmedikal.com